



## SneakPeek Early Gender DNA Test: 99.9% accuracy based on live-birth confirmation

January 12, 2021

**SneakPeek, leading provider of early gender DNA tests, ran the largest study to date of 1,029 women between 7-37 weeks into pregnancy, determining fetal sex with 99.9% accuracy.**

A large-scale clinical study of 1,029 pregnant women between 7 and 37 weeks of gestation (median = 10 weeks) was completed to determine the accuracy of the SneakPeek Early Gender DNA Test. The comparison method used was **live birth confirmation, the only true way to measure the accuracy of a test.**

Venipuncture blood samples were collected from women who were 7-37 weeks gestational age during October through December 2019. Fetal gender was determined during pregnancy using a qPCR-based SneakPeek Early Gender DNA test, and the comparison data was the actual gender of the baby at birth.

Number of Samples	1,029
Sensitivity (Male Determination)	505/506 (99.8%)
Specificity (Female Determination)	523/523 (100%)
Accuracy	1,028/1,029 (99.9%)
Gestational Age Range of Participants	7 weeks to 37 weeks (median = 10 weeks)

Of the 1,029 samples, **1,028 SneakPeek test results matched the baby's gender at birth for an overall accuracy of 99.9%.** Y-chromosome DNA was detected in 505 of the 506 samples from mothers that gave birth to baby boys, providing 99.8% sensitivity. Gender of all 523 newborn girls was identified correctly, resulting in 100% specificity. There were no cases of male DNA contamination identified in this study (i.e. no false positives).

In this study including women from 7-37 weeks into pregnancy with a median gestational age of 10 weeks, SneakPeek Early Gender DNA Test was 99.9% accurate, making it the **most dependable and reliable early fetal sex test available on the market.**

As the leading provider of early gender DNA tests since 2014, SneakPeek has been trusted by over 300,000 moms and is recommended by top physicians and obstetricians everywhere.

**SneakPeek believes in robust clinical studies that hold up under the review of impartial experts.**

That's why we submitted this data to a scientific journal for independent evaluation and publication. We're proud to announce it was published by the International Journal of Pregnancy & Child Birth on Dec. 29, 2020. See the full study here! [SneakPeek Early Gender DNA Test 99.9% Accurate](#)

